


Tri Map App Operation

RLB PI-RADS Gleason DRE	RMB PI-RADS Gleason DRE	RTZB PI-RADS Gleason	LTZB PI-RADS Gleason	LMB PI-RADS Gleason DRE	LLB PI-RADS Gleason DRE
RLM PI-RADS Gleason DRE	RMM PI-RADS Gleason DRE	RTZM PI-RADS Gleason	LTZM PI-RADS Gleason	LMM PI-RADS Gleason DRE	LLM PI-RADS Gleason DRE
RLA PI-RADS Gleason DRE	RMA PI-RADS Gleason DRE			LMA PI-RADS Gleason DRE	LLA PI-RADS Gleason DRE



The Tri Map is laid out like a twelve core biopsy
In addition it has 4 spots for transition zone sample

RLB		
PI-RADS	Gleason	DRE
RLM		
PI-RADS	Gleason	DRE
RLA		
PI-RADS	Gleason	DRE

Click on the
Right Lateral Apex Box


RLB		
PI-RADS	Gleason	DRE
RLM		
PI-RADS	Gleason	DRE
RLA		
PI-RADS	Gleason	DRE

A drop down button
appears with options

RLB		
PI-RADS	Gleason	DRE
RLM		
PI-RADS	Gleason	DRE
RLA		
PI-RADS	Gleason	DRE
PR5		


You select PI RADS 5
as MRI indicates
The box color codes
PR5 purple

RLB	PI-RADS	Gleason	DRE
RMB	PI-RADS	Gleason	DRE
RTZB	PI-RADS	Gleason	
LTZB	PI-RADS	Gleason	
LMB	PI-RADS	Gleason	DRE
LLB	PI-RADS	Gleason	DRE
RLM	PI-RADS	Gleason	DRE
RMM	PI-RADS	Gleason	DRE
RTZM	PI-RADS	Gleason	
LTZM	PI-RADS	Gleason	
LMM	PI-RADS	Gleason	DRE
LLM	PI-RADS	Gleason	DRE
RLA	PI-RADS	Gleason	DRE
PR5			
RMA	PI-RADS	Gleason	DRE
LMA	PI-RADS	Gleason	DRE
LLA	PI-RADS	Gleason	DRE



The purple scored region is now projected to the Tri Map figures

RLB	PI-RADS	Gleason	DRE
RMB	PI-RADS	Gleason	DRE
RTZB	PI-RADS	Gleason	
LTZB	PI-RADS	Gleason	
LMB	PI-RADS	Gleason	DRE
LLB	PI-RADS	Gleason	DRE
RLM	PI-RADS	Gleason	DRE
RMM	PI-RADS	Gleason	DRE
RTZM	PI-RADS	Gleason	
LTZM	PI-RADS	Gleason	
LMM	PI-RADS	Gleason	DRE
LLM	PI-RADS	Gleason	DRE
RLA	PI-RADS	Gleason	DRE
PR5	7(4+3)	T2a	
RMA	PI-RADS	Gleason	DRE
LMA	PI-RADS	Gleason	DRE
LLA	PI-RADS	Gleason	DRE



As you check and score each box the scores project to the Tri Map

RLB	PI-RADS	Gleason	DRE
	PR5	9+	T3

RMB	PI-RADS	Gleason	DRE

RTZB	PI-RADS	Gleason

LTZB	PI-RADS	Gleason

LMB	PI-RADS	Gleason	DRE

LLB	PI-RADS	Gleason	DRE

| **RLM** | PI-RADS | Gleason | DRE |
| | PR5 | 9+ | T3 |

RMM	PI-RADS	Gleason	DRE
	PR5	9+	T3

RTZM	PI-RADS	Gleason

LTZM	PI-RADS	Gleason

LMM	PI-RADS	Gleason	DRE


LLM	PI-RADS	Gleason	DRE

| **RLA** | PI-RADS | Gleason | DRE |
| | PR5 | 9+ | T3 |

RMA	PI-RADS	Gleason	DRE
	PR5	9+	T3

LMA	PI-RADS	Gleason	DRE

LLA	PI-RADS	Gleason	DRE



This is an example of a patient with extremely aggressive disease on one side. This is invaluable as you plan brachytherapy. We would make sure this side was heavily loaded to boost dose where most needed.